

PRODUCT DATA SHEET



No. CSS001

Cell Strainers

Description

Cell Strainers, Nylon mesh, Blister packed, Sterile

Purpose


Cell strainers are sterile, rapid, easy-to-use devices for isolating primary cells to consistently obtain a uniform single-cell suspension from tissues

Materials

Frame: PP (Polypropylene)

Mesh: Nylon

Dimensions (Unit: mm)







Features

- Available in 3 mesh pore sizes: 40μm, 70μm and 100μm
- · 3 different colors: blue, white and yellow, for easy identification
- · Made of a strong nylon mesh with evenly spaced mesh pores
- · Improved uniformity of single cell suspensions
- · The extended lip on the strainer enables aseptic handling with forceps
- · Design to fit perfectly into a JET 50mL centrifuge tube
- · Non-cytotoxic
- · DNase/RNase free and Non-pyrogenic
- Sterilized by irradiation SAL10⁻⁶ (ISO11137)
- . Shelf Life: 3 years after month of production
- · Manufactured in a class 100,000 room environment
- · Manufactured under ISO13485 and ISO9001 quality management system

Cell Strainers

Cat. No.	Description	Sterilization	Qty/Pk	Qty/Box	Qty/Case
CSS013040	40μm Cell Strainers, Blue, Blister packed	Yes	1	50	200
CSS013070	70µm Cell Strainers, White, Blister packed	Yes	1	50	200
CSS013100	100µm Cell Strainers, Yellow, Blister packed	Yes	1	50	200